

# Pluripotent stem cell derived organoids

## Media recipe quick reference guide

**Qkine**



PIO-070223-v1

### Three steps for choosing your growth factors

- ▶ Evidence of protein quality
  - Quantitative bioactivity data with EC50
  - Clear SDS-PAGE gel, high protein loading, good staining
  - Purity data such as mass spec, analytical reverse phase and endotoxin testing (with stringent limit <0.05. See our discussion on E.coli vs mammalian expressed growth factors)
- ▶ Good support and handling advice
- ▶ Reliable supplier, rapid delivery and complete product data

Find out more about protein quality control and growth factor selection [qkine.com/choosing-growthfactors](https://qkine.com/choosing-growthfactors)

### Quick handling guide



### Reconstitution calculator

$$\frac{\text{Mass in vial} (\mu\text{g})}{1\mu\text{g} = 1000\text{ng}} \div \frac{1\mu\text{g} = 1000\text{ng}}{\text{Desired concentration} (\mu\text{g/ml})} \times 1000 = \text{Volume to add} (\mu\text{l})$$

### How Qkine is improving growth factors for organoids

- Animal-free**  
Unmatched quality and reliability. All our proteins are made in a dedicated animal-free laboratory in Cambridge, UK.
- Total-transparency**  
Know what you're giving your cells. Stringent purity and bioactivity data for all proteins.
- Protein innovation**  
Solving stem cell culture challenges with optimized forms and animal-free firsts.

# Adult stem cell derived organoids

## Media recipe quick reference guide

**Qkine**



PAO-070223-v1.1

### Three steps for choosing your growth factors

- ▶ Evidence of protein quality
  - Quantitative bioactivity data with EC50
  - Clear SDS-PAGE gel, high protein loading, good staining
  - Purity data such as mass spec, analytical reverse phase and endotoxin testing (with stringent limit <0.05. See our discussion on E.coli vs mammalian expressed growth factors)
- ▶ Good support and handling advice
- ▶ Reliable supplier, rapid delivery and complete product data

Find out more about protein quality control and growth factor selection [qkine.com/choosing-growthfactors](http://qkine.com/choosing-growthfactors)

### Quick handling guide

- 1 Centrifuge lyophilized protein
- 2 Add reconstitution solution final concentration > 50 µg/ml
- 3 Wait 5 minutes and mix gently
- 4 Make single use aliquots and freeze -20°C or -80°C < 1 year

### Reconstitution calculator

$$\frac{\text{Mass in vial} (\mu\text{g})}{\text{Desired concentration} (\mu\text{g/ml})} \times 1000 = \text{Volume to add} (\mu\text{l})$$

1mg = 1000µg    1µg = 1000ng

### How Qkine is improving growth factors for organoids

- Animal-free  
Unmatched quality and reliability. All our proteins are made in a dedicated animal-free laboratory in Cambridge, UK.
- Total-transparency  
Know what you're giving your cells. Stringent purity and bioactivity data for all proteins.
- Protein innovation  
Solving stem cell culture challenges with optimized forms and animal-free firsts.

[qkine.com/organoids](http://qkine.com/organoids)