

- mouse - 94%
- rat - 97%
- bovine - 99%
- porcine - 99%

Product Information

- Micro-encapsulated high purity protein
- 48 mini DISC™ devices per pack
- Lot-to-lot consistency in bioactivity

Reconstitution instructions

- StemCultures

Featured applications

- Expansion of induced pluripotent, embryonic and mesenchymal stem cells

Original product page: <https://qkine.com/product/stembeads-qkine-fgf-2-154-aa-disc-device-mini/>

PDF generated: 10 May 2026

Copyright © 2026 by Qkine Ltd. All rights reserved including graphics and images.